

^{14}C Respiration/Production Methods

Mario Muscarella

22 September, 2015

Overview

The simultaneous measurement of bacterial respiration and production will give the best results for estimating growth efficiency. Here, we have constructed custom bacterial respiration/production vials (Fig 1). We will use these vials with ^{14}C -resources. However, first we must optimize the methods.




Figure 1: Bacterial Respiration and Production Vial

Optimizations

As we start there are a few things we need to optimize. We are going to start by determining the efficiency of our apparatus. We are doing this using ^{14}C -Bicarbonate. But we need to check a few things first:

1. How much Bicarbonate should we add
2. What is the efficiency of the trap

1) Bicarbonate

I did a simple experiment where I added known volumes of the bicarbonate stock to 1500 μL of scintillation cocktail (1, 5, and 10 μL). Based on this I discovered that the LSC could only measure 1 μL max. The CPM for this voucher was $\sim 2,300,000$ cpm.

2) Trap Efficiency